



CiDehol[®] ST99

Sterile 99% IPA Solution

Technical Data Sheet

- BRIEF DESCRIPTION -

CiDehol ST 99 is a sterile 99% (v/v) isopropyl alcohol solution made to clean surfaces and equipment found in cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, non-pyrogenic alcohol solution is required.

CiDehol ST 99 is filtered to 0.2 microns, bottled and double-bagged within a CLASS 100 cleanroom. It is gamma-irradiated to a sterility assurance level (SAL) of 10⁻⁶. Each lot is tested for sterility via USP 14-day sterility test and for bacterial endotoxins prior to release. A copy of lot specific documentation, which details QC, irradiation, sterility, and bacterial endotoxin limits is attached to every case of CiDehol ST 99. All lots of Cidehol ST 99 are guaranteed to contain less than 0.25 EU/mL.

- USE PROTOCOL -

Use in well-ventilated area. Hold spray bottle six to eight inches from surface and spray area until it is completely covered with solution. Or wipe CiDehol ST 99 onto surface using cleanroom grade wipers. Allow to air dry.

- INGREDIENTS -

99% Isopropyl Alcohol
1% Purified Water, USP
by Volume

- SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date can be found on the product. Do not store at temperatures above 120°F.

- PRECAUTIONARY STATEMENT -

Flammable - do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

- PACKAGING -

16 oz. Trigger Spray / 12 per case
1 Gallon / 4 per case

- CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	16 oz.	8216
	1 Gal.	8201

**For Additional Technical Information call
Decon Customer Service**

(800) 332-6647, or

See our website at www.deconlabs.com

®Decon Labs, Inc.

Revised 2/26/2014