

DECON LABORATORIES, INC.



CiDehol[®] ST (Gamma Irradiated)

For Use in Clean Rooms and
Controlled Areas
Hard Surface Disinfectant

Active ingredients:
Isopropyl Alcohol CAS# 67-63-0 70.0%
Other ingredients:
*Water 30.0%
Total: 100.0%
*USP Water for injection



**KEEP OUT OF REACH OF CHILDREN
DANGER – FLAMMABLE**

Keep away from heat or open flame.
Store at room temperature below 48.8° C (120° F)

See panel below for first aid and additional precautions

**US EPA Reg. No.: 56753-1
US EPA Est. No.: 67612-PA-001**



**NET
CONTENTS
354 mL/
12 oz.
CAT. # 8312**

**KEEP THIS INSERT
WITH CASE**

DIRECTIONS FOR USE Read the label before using

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is **NOT** to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

USE: Disinfects and sanitizes clean rooms and controlled areas such as those used in healthcare institutions, biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing facilities. Use on hard non-porous inanimate surfaces in

aseptic filling and gowning rooms, general manufacturing areas, or on machinery, tables, counters, laminar flow benches, floors, walls, carts, shelves made of hard, non-porous materials such as glass, plastic, vinyl, chrome, and stainless steel. It is compatible with most hard surface materials.

DISINFECTION: Use this product while wearing personal protection. Wear protective eyewear such as goggles or face shield, protective clothing and chemically resistant gloves. Use with adequate ventilation and spray away from eyes and face. CiDehol[®]ST is ready to use. Do not dilute. Pre-clean surface or item of heavy soil or gross filth before application. Hold container upright 6–8 inches from surface. Spray product on surface and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes wipe dry with sterilized cloth, wiper, if needed.

CiDehol[®]ST is bactericidal. The results of AOAC Efficacy Tests show that CiDehol[®]ST is an effective hard surface disinfectant against *Pseudomonas aeruginosa* (ATCC# 15442), *Staphylococcus aureus* (ATCC#6538), and *Salmonella enterica* (ATCC# 10708), in the presence of 5% fetal bovine serum organic load in 10 minutes.

SANITIZATION: Use 5 minute contact time on hard, non-porous surfaces (made from: stainless steel, glass, Corian[™], laminate flooring, aluminum, pvc, mild steel, polypropylene), in isolators, clean rooms, sterile and controlled areas, gowning rooms and on process lines, gloves, work surfaces, tools and equipment.

CiDehol[®]ST is an effective sanitizer on inanimate, non-food contact surfaces against *Staphylococcus aureus* (ATCC#6538) and *Enterobacter aerogenes* (ATCC#13048) in the presence of 5% fetal bovine serum organic load in 5 minutes.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal

Storage: Store in the original container in a dry cool place below 48.8° C (120° F).

Disposal: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control agency, or the Hazardous waste representative at the nearest EPA Regional Office for guidance.

Container Disposal: Non-refillable container. Do not reuse this container to hold materials other than pesticides or diluted pesticides (rinsate). After emptying and cleaning, it may be allowable to temporarily hold rinsate or other pesticide-related materials in the container. Contact your state regulatory agency to determine allowable practices in your state. Offer for recycling if available or puncture and dispose in a sanitary landfill, or by other procedures approved by state and local authorities. If rinsate cannot be used, follow pesticide disposal instructions.

FIRST AID

If in eyes	<ul style="list-style-type: none"> Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first five minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.
If swallowed	<ul style="list-style-type: none"> Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or Doctor. Do not give anything to an unconscious person.
If inhaled	<ul style="list-style-type: none"> Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth if possible. Call poison control center or doctor for treatment advice.
If on skin or clothing	<ul style="list-style-type: none"> Take off contaminated clothing Rinse skin immediately with plenty of water for 15 - 20 minutes. Call a poison control center or doctor for treatment advice.

Have product container or label with you when calling poison control or doctor, or going for treatment.

For emergency information on CiDehol[®]ST, call the National Pesticides Information Center at 1-800-858-7378, 6:30–4:30 pm Pacific time (PT), seven days a week. During other times, call the poison control center 1-800-222-1222.

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER: Causes irreversible eye damage. Harmful if swallowed or inhaled. Do not get in eyes or on clothing. Avoid breathing vapor or spray mist. Avoid contact with skin or clothing. Wear protective eyewear such as goggles or face shield, protective clothing and chemically resistant gloves. Wash thoroughly with soap and water after handling, before eating, drinking, chewing gum or using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

PHYSICAL AND CHEMICAL HAZARDS: FLAMMABLE liquid and vapor. Keep away from heat or open flame. Exposure to temperatures above 48.8° C (120° F) may cause bursting. Do not mix this product with other chemicals including; acids, oxidizing materials, combustible materials, halogens, peroxides, bases and metal salts.



Manufactured by:
Decon Labs, Inc.
460 Glennie Circle
King of Prussia, PA 19406