



## Technical Data Sheet

---

### **-BRIEF DESCRIPTION-**

Bathbrite is for use in water baths, incubators and other temperature controlled lab equipment. 10 to 15 drops of Bathbrite per gallon of water will keep water clear and fresh and extend the time between water replacements.

### **-PACKAGING-**

**8oz/6 per case**

### **-CATALOG NUMBER-**

<u>Item</u>	<u>Cat. No</u>
<b>6 X 8 oz.</b>	<b>3408</b>

**For Additional Technical Information call  
Decon Labs' Customer Service (800) 332-  
6647 or Visit our website at  
[www.deconlabs.com](http://www.deconlabs.com)**

Foxwest Sales.

10/09/2015