



# CiDehol<sup>®</sup>60

## Filtered 60% IPA Solution

### Technical Data Sheet

#### - BRIEF DESCRIPTION -

*CiDehol 60 is a 60% (v/v) isopropyl alcohol solution that is filtered to 0.22 microns and packaged in ready-to-use 1 gallon size. CiDehol 60 can be used for filter integrity testing or to clean surfaces in microbiology labs and tissue culture labs. It is also effective in cleaning surfaces and process equipment in pharmaceutical, biotechnology, and medical device facilities. IPA solution evaporates away completely so no rinse is required.*

#### - USE PROTOCOL -

Use in well-ventilated area. Wipe solution on surface to be cleaned until it is completely covered. Allow to air dry.

#### - INGREDIENTS -

60% Isopropyl Alcohol, USP  
40% Purified Water, USP

#### - SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date can be found on the product. Do not store at temperatures above 120°F.

#### - PRECAUTIONARY STATEMENT -

Do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

#### - PACKAGING -

1 Gallon / 4 per case

#### - CATALOG NUMBER -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	1 Gal.	8800

**For Additional Technical Information call  
Decon Customer Service  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)**

©Decon Labs, Inc.

Revised 10/24/2013