



CiDehol[®] 70 Wipes

Technical Data Sheet

- BRIEF DESCRIPTION -

CiDehol 70 wipes are cleanroom wipes that are nonwoven polyester/cellulose wipes pre-saturated with 70% isopropyl alcohol. Each pack is in a convenient, canister containing 160 durable wipes, 6 x 7 inches.

- DIRECTIONS FOR DISPENSING -

The C-folded wipers are made of a polyester cellulose, unwoven and are easily drawn from canister and limits use to one wiper at a time.

- USE PROTOCOL -

These wipes are designed for use in cleanroom areas within pharmaceutical, biotechnology and medical device facilities.

- ACTIVE INGREDIENTS -

Isopropyl Alcohol, USP Grade (CAS #67-63-0)....70 %
Deionized Water (CAS #7732-18-5)...30 %

- STORAGE AND DISPOSAL -

Store in cool area. Do not reuse wipes. Dispose of used wipe and empty dispenser in trash.

- PRECAUTIONARY STATEMENT -

CAUTION: FLAMMABLE. Do not expose wipes to heat, open flame, sparks, or elevated temperatures.

- FIRST AID -

For all cases of exposure consult with a Physician.
Eye Contact: Immediately flush with water for at least 15 minutes.
Skin Contact: Immediately rinse with water.
Inhalation: Immediately remove person to fresh air.
Ingestion: Immediately provide water to dilute ingested material.

- PACKAGING -

160 wipes per canister / 12 canisters per case

- CATALOG NUMBER -

| <u>Manufactured for:</u> | <u>Item</u> | <u>Cat. No.</u> |
|--------------------------|-------------|-----------------|
| Decon Labs, Inc. | 160 Wipes | 8775 |

**For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

®Decon Labs, Inc.

Revised 08/18/2011