

DLI

KOPECTE PURE ETHANOL 190 PROOF



KOPECTE 190 proof ethanol (95% ethanol) is a pure ethanol product made to USP specifications and used for research, healthcare, industrial, or production applications as a cleaner, solvent or reagent.

Additional applications include food processing: Decon Labs Inc.'s 190 Proof pure Ethanol meets the standards listed in the food chemicals codex (FCC) 6th edition ethyl alcohol monograph and is considered "food grade."

It is available in convenient packaging configurations from stock. Custom solutions or packaging configuration may be available.

Please note that pure ethanol is regulated by the US ATF/TTB and depending on the taxable status of your institution, pure ethanol products are taxed according to the government schedule of taxes detailed in 27 CFR Chapter 1 section 27.41.

To purchase tax-free - you must have an Industrial Alcohol User Permit from the AFT/TTB, specifically marked for "the withdraw and use of alcohol free of tax" Tax free pure ethanol is intended for industrial use only. To apply for this permit visit www.ttb.gov and search for "Form 5150.22" - which is the application for the Industrial Alcohol User Permit.

INGREDIENTS

95% Ethyl Alcohol (CAS #64-17-5)

5% Purified Water (CAS #7732-18-5)

STORAGE: Do not store Ethanol at temperatures above 120° F.

Catalog Numbers / Packaging

Item	CASE QTY	Decon # Tax-Free	Decon # Tax-Paid
190 Proof, 1 Pint	24	V1116	V1116TP
190 Proof, 1 Gallon (plastic)	4	V1101	V1101TP
190 Proof, 1 Gallon (glass)	4	V1101G	V1101GTP
190 Proof, 5 Gallon (plastic w/ spout)	-	V1105	V1105TP
190 Proof, 5 Gallon (Hedwin Cubitainer w/ spigot)	-	V1105HC	V1105HCTP
190 Proof, 5 Gallon (metal w/ spout)	-	V1105M	V1105MTP
190 Proof, 5 Gallon (w/ a spigot closure*)	-	V1105SG	V1105SGTP
190 Proof, 55 Gallon (plastic drum)	-	V1155	V1155TP
190 Proof, 55 Gallon (metal drum)	-	V1155M	V1155MTP

*Spigot available separately

If you possess an Industrial Alcohol Users Permit order using the Decon Tax-Free part number. All others would order using the Tax-Paid part number.

Rev. 12/2022