



**200 Proof Ethanol  
100 % Ethanol**

**Technical Data Sheet**

**- BRIEF DESCRIPTION -**

KOPTEC 200 proof ethanol (100% ethanol) is a pure ethanol product made to USP specifications and used for research, healthcare, industrial, or production applications as a cleaner, solvent, or reagent.

It is available in convenient packaging configurations from stock. Custom solutions or packaging configuration may be available.

Please note that pure ethanol is regulated by the US ATF/TTB and depending on the taxable status of your institution, pure ethanol products may be taxed according to the government schedule of taxes detailed in 27 CFR Chapter 1 section 27.41.

To purchase Ethanol free of Federal Excise Tax,(tax free), you must have an Industrial Alcohol User Permit from the AFT/TTB, specifically marked for “the withdraw and use of alcohol free of tax” Tax free pure ethanol is intended for industrial use only. To apply for this permit visit [www.ttb.gov](http://www.ttb.gov) click Laboratory, click Laboratory Forms- click “Form 5150.22”

**- INGREDIENTS-**

100 % Ethyl alcohol      CAS Number 64-17-5

**- STORAGE -**

Do not store Ethanol at temperatures above 120 degrees F.

**- PACKAGING -**

**1 Pint or 1 Glass Pint / 24 per case  
1 Gallon or 1Glass Gallon / 4 per case  
5 Gallons / each  
55 Gallons / each**

**- CATALOG NUMBERS -**

<u>Manufacturer</u>	DLI Labs, Inc.
<u>Item</u>	<u>Cat. No.</u>
24x1 Pint	V1016
24x1 Pint, Glass	V1016G
4x1 Gallon	V1001
4x1 Gallon, Glass	V1001G
5 Gallon	V1005
Gallon Hedwin Cubitainer	V1005HC
5 Gallon Metal Container	V1005M
5 Gallon Spigot Closure	V1005SG
55 Gallon Drum	V1055
55 Gallon Metal Drum	V1055M

*For tax-paid items add “TP” to the end of the catalog number. For example, the item tax-paid 24 X 1 pints has a catalog number V1016TP.*

**For Additional Technical Information call  
Customer Service (800) 332-6647**

*Revised -5/19/2016*