



LeStab®

Inoculating Needles

Technical Data Sheet

- BRIEF DESCRIPTION -

Decon inoculating needles are made of high quality, nichrome wire which is designed to cool rapidly after flaming. Most of our needles are permanently attached to lightweight, insulated handles which are perfectly balanced for easy manipulation. The nichrome wire cools rapidly after flaming.

Decon also offers 3 inch inoculating needles. These needles fit in Decon's LeLoop Brass holders. The Brass Holders come in 6 inch and 8 inch size.

- USE PROTOCOL -

LeStab inoculating needles with attached handles can be used and flamed many times before degrading the nichrome wire. Once degradation occurs, dispose of the needle with the attached handle and use a new one.

When the 3 inch inoculating needle needs to be replaced, remove the needle from the brass handle. Dispose of the used needle and place a fresh needle on the Brass holder.

- CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs	5 pack Silver Handle attached 11 inch w/handle	MP186-5
	1 Each Blue Angled Handle attached 8 inch w/handle	MP 4010-1G
	1 Each White Handle attached 8 inch w/handle	MP 4000-1G
	12 pack 3 inch inoculating needles only Brass Holders sold separately	MP 3325-12G
	10 pack Le Loop Brass Holder 6 inch 8 inch	MP6-10 MP8-10

**For Additional Technical Information call
Decon Customer Service Hotline
(800) 332-6647, or
See our website at www.deconlabs.com**

©Decon Labs, Inc.

Revised 02/19/2016