



## Technical Data Sheet

### -BRIEF DESCRIPTION -

*Peroxygen ST is a sterile 6% (v/v) Hydrogen Peroxide solution made with Water for Injection (WFI). This ready-to-use solution can be used to clean surfaces and equipment found in cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, hydrogen peroxide solution is required.*

*Peroxygen ST is filtered to 0.2 microns, bottled and double-bagged within a CLASS 100 cleanroom. It is gamma-irradiated to a sterility assurance level (SAL) of 10<sup>-6</sup>. Each lot is tested for sterility via USP 14-day sterility test. A copy of lot specific documentation, which details QC, irradiation, and sterility, is attached to every case of Peroxygen ST.*

### - INGREDIENTS -

6% (v/v) Hydrogen Peroxide  
94% (v/v) USP Water for Injection

### - SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date is on the product.

### - PRECAUTIONARY STATEMENT -

Causes serious eye irritation. Wear eye protection/face protection

### - PACKAGING -

16 oz. (473mL) Trigger Spray / 12 per case

### - CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	16 oz. Bottle	9716

For Additional Technical Information call  
Decon Customer Service  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)

©Decon Labs, Inc.

12/21/2016