

## Technical Data Sheet

---

### - BRIEF DESCRIPTION -

*SprayKlean™ is a ready-to-use all-purpose cleaner with a fresh citrus scent. SprayKlean™ can be used on any washable surface and is packaged in a convenient 32oz. trigger spray bottle.*

### - APPLICATIONS -

While SprayKlean™ is safe to use on most surfaces, always test the product on an unexposed area of surface to be treated to ensure desired results. Good for cleaning:

- Labware
- Instruments
- Process Equipment
- Production mixers, fillers, storage containers
- Optics

### - WORKING SOLUTIONS -

SprayKlean™ is ready-to-use solution. Simply spray the surface to be treated and immediately wipe clean with a dry, lint-free cloth. Tougher stains may require a contact time of up to 3 minutes. Do not leave solution on surface for more than the recommended contact time.

### - STORAGE -

Keep contents in a clean and dry environment.

### - PACKAGING -

32oz. / 6 per case

### - CATALOG NUMBERS-

| <u>Item</u>   | <u>Cat. No.</u> |
|---------------|-----------------|
| 32oz. Trigger | 4732            |

---

**For Additional Technical Information call  
Decon Customer Service  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)**

---

®Decon Labs, Inc.      Created 08/01/2014