



Memorandum

July 1, 2022

From: Karen Wend
610-755-0800
kwend@deconlabs.com

Subject: Ethanol, Pure, 200,190,140 Proof expiration

Decon's Ethanol, 200,190, and 140 proof is a pure ethanol product made to USP specifications and used for research, healthcare, industrial, or production applications as a cleaner, solvent, or reagent.

Due to the stability of the product there is no expiration date. As a non-binding recommendation if an expiration date is needed, Decon recommends an expiration date of 3 years after the manufacturing date. However, using this product past that should not affect its performance.

Please feel free to contact me at Decon Labs if you require additional information.

A handwritten signature in blue ink that reads 'Karen Wend'.

Manager, Quality and Technical Services