



Memorandum

From: Charlotte Veloski
Manager, Quality and Technical Services
Decon Labs
460 Glennie Circle
King of Prussia, PA 19406
610-755-0800 Ext. 31
cveloski@deconlabs.com

Date: January 28, 2016
Subject: Ethanol, Pure, 200,190,140 Proof expiration

Decon's Ethanol, 200,190, and 140 proof is a pure ethanol product made to USP specifications and used for research, healthcare, industrial, or production applications as a cleaner, solvent, or reagent.

Due to the stability of the product there is no expiration date. As a non-binding recommendation if an expiration date is needed, Decon recommends an expiration date of 3 years after the manufacturing date. However, using this product past that should not affect its performance.

Please feel free to contact me at Decon Labs if you require additional information.

Charlotte Veloski