

Alcohol
Hydrogen Peroxide &
Purified Water

70% Medical Device
6% Hydrogen Peroxide
Research
Process Equipment
Biotechnology
Wipes
IPA
Healthcare
Laboratory
60% EtOH
Cleanroom
70% EtOH
Ready-to-use
Cleanroom
Alcohol
IPA
Laboratory
Medical Device
Research
Process Equipment
Water
Pharmaceutical
Process Equipment
Ready-to-use
Research
IPA 60%
6% Hydrogen Peroxide



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com

ALCOHOL

CiDehol® 70

70% Isopropyl Alcohol

CiDehol® 70 is a ready-to-use IPA solution that is filtered to 0.2 µm and packaged in convenient trigger spray and gallon bottles. Ideal for cleaning lab surfaces, production areas, and equipment. Evaporates rapidly – no rinsing, no residues.



CiDehol® 60

60% Isopropyl Alcohol

CiDehol® 60 is a ready-to-use IPA solution that is filtered to 0.2 µm and packaged in convenient gallon bottles. Ideal for filter integrity testing or for cleaning lab surfaces, production areas, and equipment. Evaporates rapidly – no rinsing, no residues.



Sanihol® 70

70% Denatured Ethyl Alcohol

Sanihol® 70 is a ready-to-use 70% solution of denatured ethanol for cleaning surfaces in labs and production areas. Perfect for use in tissue culture hoods and microbiology labs – it evaporates away leaving no residues. Filtered to 0.2 µm to remove particulate contamination.



CiDehol® 70 Wipes

Pre-Saturated IPA Wipes

CiDehol® 70 Wipes are pre-saturated with 70% isopropyl alcohol and 30% USP grade water. These wipes are ideal for cleaning hard surfaces, production areas and equipment in labs, cleanrooms and other facilities. Simply wipe the surface to be cleaned and allow excess moisture to evaporate. The wipes come 175 to a canister.



PURIFIED WATER

AquaPur™ ST

Purified Water USP Grade

AquaPur™ ST is sterile purified water for use in dilution of sterile disinfectants, manufacturing, general laboratory and cleanroom applications. Purified water is filtered to 0.2 µm, irradiated to a SAL 10⁻⁶ and tested for sterility and endotoxin using current USP method. Meets USP specs. Each case is shipped with a Lot Specific Document. Not for parenteral use. Not for injection, irrigation or infusion.



AquaPur™

Purified Water Meets USP/EP Specs

AquaPur™ Purified Water has been tested to meet USP/EP specifications and packaged in convenient point of use packaging, for use in laboratory or production applications requiring Purified Water. The product is filtered to 0.2µm. Guaranteed to meet specifications through expiration date.



PRODUCT GUIDE

PRODUCT	SIZE	DECON CAT. #	CASE SIZE
CiDehol® 70 70% IPA	16 oz. Trigger Spray	8416	12 x 16 oz.
	32 oz. Trigger Spray	8432	12 x 32 oz.
	1 gal.	8401	4 x 1 gal.
	5 gal.	8406	NA
	55 gal.	8405	NA
CiDehol® 60 60% IPA	1 gal.	8800	4 x 1 gal.
Sanihol® 70 70% Denatured Ethanol	16 oz. Trigger Spray	8616	12 x 16 oz.
	1 gal.	8601	4 x 1 gal.
CiDehol® 70 Wipes IPA Surface Wipes	175 Wipes	8776	12 x 175 Wipes
AquaPur™ ST Purified Water USP Grade	1 gal.	8904	4 x 1 gal.
AquaPur™ Purified Water Meets USP/EP Specs	1 gal.	9301	4 x 1 gal.
	5 gal.	9305	EACH
	55 gal.	9355	EACH