

# Detergents & Cleaners

**Acidic**  
**Powder**  
**Liquid**

**Ultrasonics**

**Neutral**  
Water Bath Conditioner  
**Liquid**

**Enzymatic** Ready-to-use  
**Detergents**

**Acidic**  
**Enzymatic**

**Cleaners** **Ultrasonics**

**Alkaline**  
Enzymatic Ready-to-use  
**Lens Cleaner**

**Detergents**

Concentrated **Manual**  
**Cleaners**  
Concentrated  
**Ready-to-use**

**Water Bath Conditioner**  
**Lens Cleaner**  
**Concentrated**



460 Glennie Circle  
King of Prussia, PA 19406  
(800) 332-6647  
[www.deconlabs.com](http://www.deconlabs.com)

## MANUAL CLEANERS

### Contrad® 70

#### *The Ultimate Labware Cleaner*

Contrad® 70 is a liquid detergent concentrate that gets labware clean by soaking, scrubbing and ultrasonics. It rinses away totally, leaves zero residues. Phosphate-free, disposable straight to drain, Contrad® 70 is perfect for all critical cleaning, from labware to process equipment, ultrasonic cleaning and radioactive decontamination. The original and still the best... use it and see!



### Dri-Clean®

#### *Powdered Detergent*

Dri-Clean® is a powdered detergent for soak, manual, and ultrasonic cleaning applications. User friendly bottle provides easy dispensing and eliminates caking. Completely soluble, free-rinsing, and residue-free. Versatile formula removes a wide variety of tenacious contaminants.



### Neutrad®

#### *Neutral pH Liquid Cleaner*

Neutrad® is a liquid detergent concentrate for hand scrubbing your labware which offers high-performance cleaning with the safety of a near-neutral pH. Phosphate-free, totally rinsable, disposable straight to drain, Neutrad® is the perfect replacement for old-fashioned caustic soap powders.



### Citrad®

#### *Acidic Liquid Detergent*

Citrad® is a citric acid-based detergent concentrate ideal for removing salts, scale, and metal-oxide. Also used for metal brightening and pharmaceutical-process equipment. Completely phosphate-free, totally free-rinsing and biodegradable.



### Contrex® AL

#### *Alkaline Liquid Detergent*

Contrex® AL is a concentrated alkaline liquid detergent formulated for use in manual, soak or ultrasonic cleaning of labware, scientific apparatus and process equipment. Contrex® AL is free-rinsing, phosphate-free and disposable straight to drain.



### Contrex® AP

#### *Alkaline Powdered Detergent*

Contrex® AP is a concentrated powdered labware detergent formulated for use in manual, soak or ultrasonic cleaning of labware, scientific apparatus and process equipment. Contrex® AP is completely soluble, free-rinsing and safe to dispose straight to drain. The user friendly bottle allows for easy dispensing and eliminates caking.



### DeSCAL®

#### *Acidic Neutralizing Rinse (Liquid Descaler)*

DeSCAL® is an inorganic acid-based detergent and descaler concentrate for removing hard water mineral deposits, rust and crusty build up from labware, water baths and stills. Low-foaming formula allows it to be used for automatic washing applications.



### Contrex® CA

#### *Acidic Liquid Detergent*

Contrex® CA is a concentrated acidic liquid detergent made for use in descaling, mineral removal, oxide removal or in applications where an acidic detergent is indicated. Formulated with organic acid, Contrex® CA can be used as a manual, soak or ultrasonic cleaner.



## AUTO WASHERS

### Dri-Contrad®

#### *Powdered Autowasher Detergent*

Dri-Contrad® is a powdered detergent that supercharges your automatic labware washer. Non-foaming, instantly soluble, rinsable, and chlorinated for biological control. Dri-Contrad is the only powder you need for your labware washer.



### Contrad® NF

#### *Liquid Non-Foaming Autowasher Detergent*

Contrad® NF is a non-foaming super-concentrated liquid detergent formulated for automatic labware washers. Phosphate-free, and totally free-rinsing, Contrad NF maximizes the efficiency of your labware machine



### Contrex® AW

#### *Powdered Autowasher Detergent*

Contrex® AW is a concentrated low-foaming powdered cleanser formulated for use in all automatic labware washing machines. Contrex AW is completely soluble, free-rinsing and biodegradable.



### Contrad® 100

#### *Liquid CIP Cleaner*

Contrad® 100 is a low-foaming detergent designed for clean-in-place systems and other high-pressure spray applications, including automatic washers. Alkaline formulation is phosphate-free, residue-free, and biodegradable.



See last page for additional information on Decon detergents.

## SPECIALTY CLEANERS

### ElimiNase®

#### RNase/DNase Remover

ElimiNase® is a non-toxic, ready-to-use solution for instantly removing RNase, DNase, unwanted DNA, and radioactive contamination from labware and equipment. Safe and effective. Eliminates need for autoclaving or using DEPC.



### Optical Lens Cleaner

#### Optical Glass Surface Cleaner

Optical Lens Cleaner is excellent for removing oil and smudges from objectives or any optical surface. It gently cleans and will not affect a coated optic or soften the mounting that holds the objective in place. This lens cleaner is laboratory-safe and non-flammable. Regular cleaning helps extend the life



### Contrex® EZ

#### Powdered Enzymatic Detergent

Contrex® EZ is a powdered detergent formulated for use as an enzymatic detergent for manual or soak cleaning of labware and equipment contaminated with proteins and fats. Contrex® EZ prevents corrosion, is free-rinsing and disposable straight to drain.



### Jet Stream®

#### High Pressure Dust Remover

Jet Stream® High Pressure Dust Remover is moisture-free propelled air used to remove microscopic dust from negatives, microscope slides, lenses, balances and other laboratory equipment.



### Contrex® CF

#### Cation / Phosphate-Free Liquid Detergent

Contrex® CF is a concentrated liquid detergent formulated for use in cleaning of equipment and electronics that are sensitive to metal cations and phosphates. Contrex® CF is free-rinsing, phosphate-free and disposable straight to drain. Can be used for soak, scrub, sonication and auto washers.



### BathBrite™

#### Water Conditioner for Water Baths

Bath Brite™ can be used to condition the water in water baths, incubators, and other temperature controlled lab equipment. 10 to 15 drops per gallon of water will keep water clear and fresh. Using Bath Brite™ will extend the time between water replacement.



## Detergents & Cleaners additional information



= Liquid



= Powder

DeSCAL®



Contrex® CA



Citrad®



Neutrad®



Contrex® AL



Dri-Clean®



Contrex® AP



Contrad® 70



Dri-Contrad®



Contrad® NF



Contrex® AW



Contrad® 100



pH: ACIDIC

NEUTRAL

ALKALINE

SOAK, MANUAL SCRUB, ULTRASONIC CLEANING

ALKALINE

AUTO WASHERS

## PRODUCT GUIDE

PRODUCT	SIZE	DECON CAT. #	CASE SIZE
<b>Contrad 70®</b> <i>Ultimate Labware Liquid Detergent</i>	1 Liter	1002	12 x 1 L
	5 Liters	1003	4 x 5 L
	20 Liters	1004	–
	210 Liters	1005	–
<b>Dri-Clean®</b> <i>Powdered Detergent</i>	2 kG	2502	6 x 2 kG
	50 kG	2501	–
<b>Neutrad®</b> <i>Neutral pH Liquid Cleaner</i>	950 mL	3003	12 x 950 mL
	5 Liters	3001	4 x 5 L
	20 Liters	3002	–
<b>Citrad®</b> <i>Acidic Liquid Detergent</i>	3.8 Liters	4401	4 x 3.8 Liters
<b>Contrex® AL</b> <i>Alkaline Liquid Detergent</i>	1 qt.	5332	12 x 1 qt.
	3.8 Liters	5301	4 x 3.8 Liters
<b>Contrex® AP</b> <i>Alkaline Powdered Detergent</i>	4 lbs.	5204	9 x 4 lbs.
	25 lbs.	5225	–
<b>DeSCAL®</b> <i>Acidic Neutralizing Rinse (Liquid Descaler)</i>	1 Liter	4001	12 x 1 L
	3.8 Liters	4004	4 x 3.8 Liters
	115 Liters	4030	–
<b>Contrex® CA</b> <i>Acidic Liquid Detergent</i>	3.8 Liters	5701	4 x 3.8 Liters
<b>Dri-Contrad®</b> <i>Powdered Autowasher Detergent</i>	2 kG	2002	6 x 2 kG
	5 kG	2005	4 x 5 kG
	12.5 kG	2007	–
	50 kG	2004	–
<b>Contrad® NF</b> <i>Liquid Non-foaming Autowasher Detergent</i>	3.8 Liters	6002	4 x 3.8 Liters
	10 Liters	6001	2 x 10 L
	115 Liters	6003	–
	210 Liters	6055	–
<b>Contrex® AW</b> <i>Powdered Autowasher Detergent</i>	4 lbs.	5604	9 x 4 lbs.
	25 lbs.	5625	–
	50 lbs.	5650	–
<b>Contrad® 100</b> <i>Liquid Clean-in-place Detergent</i>	3.8 Liters	1504	4 x 3.8 Liters
	20 Liters	1503	–
	115 Liters	1506	–
	210 Liters	1505	–
<b>ELIMINase®</b> <i>RNase/DNase Remover</i>	250 mL	1101	12 x 250 mL
	950 mL	1102	12 x 950 mL
<b>Contrex® EZ</b> <i>Powdered Enzymatic Detergent</i>	4 lbs.	5404	9 x 4 lbs.
	25 lbs.	5425	–
<b>Contrex® CF</b> <i>Cation / Phosphate-Free Liquid Detergent</i>	3.8 Liters	5801	4 x 3.8 Liters
<b>Jet Stream®</b> <i>High Pressure Dust Remover</i>	11 oz.	3711	6 x 11 oz.
<b>Optical lens cleaner™</b> <i>Optical Glass Surface Cleaner</i>	6 x 2 oz.	1100	8 x 6/pk
<b>Bath Brite™</b> <i>Water Conditioner for Water baths</i>	8 oz.	3408	6 x 8 oz.