

**Sterile Alcohol  
Sterile Disinfectants  
Sterile Water**

**70%** **Medical Device** **Research** **Phenolic** **Alcohol** **70%** **EtoH** **Sterile Water**  
**Process Equipment** **Low pH** **Concentrated**  
**Biotechnology** **Disinfectants** **Broad-Spectrum**  
**Wipes** **Disinfectants** **Alcohol** **6% Hydrogen Peroxide** **Alcohol**  
**70%** **Healthcare** **Laboratory** **Sprays** **Pharmaceutical** **70%**  
**Pathogens** **Wipes** **Wipes** **Process Equipment** **IPA** **IPA**  
**EtoH** **Cleanroom** **Ready-to-use** **EtoH**  
**IPA Water**  
**Purified Water**  
**6% Hydrogen Peroxide**



460 Glennie Circle  
King of Prussia, PA 19406  
(800) 332-6647  
[www.deconlabs.com](http://www.deconlabs.com)

## STERILE ALCOHOL

### CiDehol® ST

#### *Sterile 70% Isopropyl Alcohol*

CiDehol® ST is a ready-to-use 70% IPA solution in non-aerosol trigger spray and one gallon bottles. It is filtered to 0.2 µm, individually double-bagged and gamma-irradiated. Each lot is tested for bacterial endotoxin and sterility. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing. Made with WFI. EPA registered.



### Sanihol® ST

#### *Sterile 70% Denatured Ethyl Alcohol*

Sanihol® ST is a ready-to-use 70% Ethanol solution in non-aerosol trigger spray and one gallon bottles. It is filtered to 0.2 µm, individually double-bagged and gamma-irradiated. Each lot is tested for bacterial endotoxin and sterility. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing. Made with WFI. EPA registered.



### CiDehol® ST 99

#### *Sterile 99% Isopropyl Alcohol*

CiDehol® ST 99 is a ready-to-use 99% IPA solution in non-aerosol trigger spray and one gallon bottles. It is filtered to 0.2 µm, individually double-bagged and gamma-irradiated. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing.



## STERILE DISINFECTANTS & CLEANERS

### Peroxigen™ ST

#### *6% Sterile Hydrogen Peroxide*

Sterile, double-bagged, gamma-irradiated 6% Hydrogen Peroxide solution made with Water-For-Injection (WFI) for use in cleanrooms and critical areas requiring a sterile cleaner. No residues, low-toxicity, non-flammable, and zero VOC profile help to assist in VOC compliance. Each trigger spray bottle is double-bagged with a sterile spray head and each case label contains a QR Code that is linked to a Lot Specific Document detailing parameters of QC, irradiation, and sterility testing. Not EPA registered.



## STERILE WATER




### AquaPur™ ST

#### *Purified Water USP Grade*

AquaPur™ ST is sterile purified water for use in dilution of sterile disinfectants, manufacturing, general laboratory and cleanroom applications. Purified water is filtered to 0.2 µm, irradiated to a SAL 10<sup>-6</sup> and tested for sterility and endotoxin using current USP method. Meets USP specs. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing. Not for injection, irrigation or infusion.



## PRODUCT GUIDE

PRODUCT		SIZE	DECON CAT. #	CASE SIZE
<b>CiDehol® ST 70%</b> <i>Sterile 70% IPA</i>		12 oz. Trigger Spray	8312	16 x 12 oz.
		16 oz. Trigger Spray	8316	12 x 16 oz.
		32 oz. Trigger Spray	8332	12 x 32 oz.
		1 gal. Refill	8301	4 x 1gal.
<b>Sanihol® ST 70%</b> <i>Sterile Denatured Ethanol</i>		16 oz. Trigger Spray	8116	12 x 16 oz.
		1 gal.	8101	4 x 1 gal.
<b>CiDehol® ST 99%</b> <i>Sterile IPA</i>		16 oz. Trigger Spray	8216	12 x 16 oz.
		1 gal.	8201	4 x 1 gal.
<b>Peroxigen™ ST</b> <i>Sterile 6% Hydrogen Peroxide</i>		16 oz. Trigger Spray	9716	12 x 16 oz.
<b>AquaPur™ ST</b> <i>Sterile Water</i>		1 gal.	8904	4 x 1 gal.