



CiDehol[®] 60

60% IPA Solution



CiDehol[®] 60 is a 60% (v/v) isopropyl alcohol solution that is filtered to 0.22 microns and packaged in ready-to-use 1 gallon size. CiDehol 60 can be used for filter integrity testing or to clean surfaces in microbiology labs and tissue culture labs. It is also effective in cleaning surfaces and process equipment in pharmaceutical, biotechnology, and medical device facilities. IPA solution evaporates away completely so no rinse is required.

Use Protocol: Use in well-ventilated area. Spray or wipe solution on surface to be cleaned until it is completely covered. Allow to air dry.

Storage: Do not store at temperatures above 120°F.

Precautionary Statements

Flammable - do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

INGREDIENTS

60 % Isopropyl Alcohol, USP

40% Purified Water, USP

SHELF LIFE
2 YEARS

Size	Case Qty	Decon Cat #
1 Gallon	4	8800
5 Gallon	-	8805

Rev. 09/2022