



CiDehol 70 Wipes 70% IPA Solution



CiDehol® 70 Wipes are strong, tear-resistant, lint-free wipes made for use on hard surfaces in cleanrooms, labs and other areas where a pre-saturated cleaning wipe is a convenience. Wipes are 100% PET non-woven and are 5" x 8" (175 per canister). They come in an easy-to-use dispensing canister. These canisters have a unique slam-shut lid that helps keep wipes saturated.

Use Protocol: These wipes are designed for use in cleanroom areas with in pharmaceutical, biotechnology, and medical device facilities.

Storage & Disposal: Store in a cool area. Do not reuse wipes. Dispose of used wipes and the empty dispenser in trash.

Precautionary Statements

Flammable - do not expose wipes to heat, open flame, sparks, or elevated temperatures.

INGREDIENTS

70 % Isopropyl Alcohol, USP

30% Purified Water, USP

**SHELF LIFE
2 YEARS**

CiDehol 70 Wipes

**Wipes: 5" x 8"
175 wipes per cannister
12 Cannisters per case**

**Decon Catalog#
8776**

Rev. 11/2022