



CitriSolv®

SOLVENT & CLEARING AGENT



CitriSolv® is a d-limonene based solvent and clearing agent that can be used as a safe alternative to xylene and ethyl acetate. It is soluble with absolute alcohol and will dissolve paraffin, glues, adhesives, and mounting media. It will not harden or shrink tissue.

CitriSolv® is biodegradable and non-corrosive. It does not contain benzene and has low levels of toxicity.

USE PROTOCOL: Use for manual or automatic tissue processor staining procedures in histology/cytology applications that require

STORAGE & DISPOSOL: Store in original container with lid tightly sealed. Follow federal, state and local restrictions for disposal.

PRECAUTIONARY STATEMENTS: Keep out of reach of children.

Avoid contact with eyes, skin, and clothing. Avoid ingestion and inhalation. Empty containers retain product residue and can be dangerous. Avoid contact with heat, sparks, or flame.

FIRST AID:

Eye Contact: Immediately flush eyes with plenty of water for at least fifteen minutes, occasionally lifting the upper and lower eyelids. Contact a physician.

Skin Contact: Flush skin with plenty of soap and water for at least fifteen minutes while removing contaminated clothing and shoes. Contact a physician if irritation develops or persists.

Ingestion: Do NOT induce vomiting. If victim is conscious and alert, give 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Contact a physician immediately.

Inhalation: Remove from exposure to fresh air immediately. Contact a physician.

PACKAGING

4 X 1 GALLON

Decon Catalog No. 1601

Rev. 12/2022