



Bluing Tablets



Decon's *Bluing Tablets* are to be used to color the water in a toilet, via a blue dye, in order to prevent the water from being used to dilute a urine sample for drug testing. Packaged as two tablets per each foil packet, bluing tablets are an efficient and inexpensive means of preventing critical drug sample tampering.

DIRECTIONS FOR USE: Place one or two tablets in a toilet bowl. The effervescent action will immediately start to color the water blue.

STORAGE: Store the packets in a cool, dry place.

DO NOT INGEST. FOR ADDING DYE TO TOILET WATER ONLY.

**2 Tablets per pack
50 packs per case
Decon Catalog No. 1220**

Rev. 08/2022