



DECON'S PURE ETHANOL 140 PROOF

Decon's 140 proof ethanol (70% ethanol) is a pure ethanol product made for research, healthcare, industrial, or production applications as a cleaner, solvent, or reagent. It is available in convenient packaging configurations. Custom solutions or packaging configuration may be available.

Please note that pure ethanol is regulated by the US ATF/TTB and depending on the taxable status of your institution, pure ethanol products are taxed according to the government schedule of taxes detailed in 27 CFR Chapter 1 section 27.41.

To purchase tax-free - you must have an Industrial Alcohol User Permit from the AFT/TTB, specifically marked for "the withdraw and use of alcohol free of tax" Tax free pure ethanol is intended for industrial use only. To apply for this permit visit www.ttb.gov and search for "Form 5150.22" - which is the application for the Industrial Alcohol User Permit.

INGREDIENTS
70% Ethyl Alcohol (CAS #64-17-5)
30% Purified Water (CAS #7732-18-5)

STORAGE: Do not store Ethanol at temperatures above 120° F.

Catalog Numbers / Packaging			
Item	CASE QTY	Decon # Tax-Free	Decon # Tax-Paid
140 Proof, 32oz. Plastic Bottle	12	2432	2432TP
140 Proof, 1 Gallon Plastic Bottle	4	2401	2401TP
140 Proof, 5 Gallon Plastic Pail w/ spout	-	2405	2405TP
140 Proof, 5 Gallon Metal Pail	-	2405M	2405MTP
140 Proof, 55 Gallon Metal Drum	-	2455M	2455MTP

If you possess an Industrial Alcohol Users Permit order using the Decon Tax-Free part number. All others would order using the Tax-Paid part number.

Rev. 12/2022