



140 Proof Ethanol

70% Ethanol

Technical Data Sheet

- BRIEF DESCRIPTION -

Decon's 140 proof ethanol (70% ethanol) is a pure ethanol product made for research, healthcare, industrial, or production applications as a cleaner, solvent, or reagent.

It is available in convenient packaging configurations. Custom solutions or packaging configuration may be available.

Please note that pure ethanol is regulated by the US ATF/TTB and depending on the taxable status of your institution, pure ethanol products are taxed according to the government schedule of taxes detailed in 27 CFR Chapter 1 section 27.41.

To purchase you must have an Industrial Alcohol User Permit from the AFT/TTB, specifically marked for "the withdraw and use of alcohol free of tax" Tax free pure ethanol is intended for industrial use only. To apply for this permit visit www.ttb.gov click Laboratory, click Laboratory Forms-click "Form 5150.22"

- INGREDIENTS-

70 % Ethyl alcohol CAS # 64-17-5
30 % Purified Water CAS # 7732-18-5

- STORAGE -

Do not store Ethanol at temperatures above 120 degrees F.

- PACKAGING -

32 oz/ 12 per case
1 Gallon / 4 per case
5 Gallons/each
55 Gallons/each

- CATALOG NUMBERS -

<u>Manufacturer</u>	Decon Labs, Inc.
<u>Item</u>	<u>Cat. No.</u>
32 oz	2432
1 Gallon(plastic)	2401
5 Gallon (plastic)	2405
5 Gallon Metal Container	2405M
55 Gallon Metal Container	2455M

For tax-paid items add "TP" to the end of the catalog number. For example, the item tax-paid 4 X 1 gallon has a catalog number 2401TP.

**For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

Decon Labs, Inc. Revised 5/19/2016