



**Optical
Lens Cleaner**



*Decon's **Optical Lens Cleaner** is a gentle cleaning solution that can be used to clean all optical glass surfaces including microscopes, optics, adsorption cells, and lenses.*

***Optical Lens Cleaner** contains no abrasives or silicones so it will not etch, scratch, streak, or leave interfering residues. Safe to use with cemented lenses.*

APPLICATIONS: Apply to optical glass surface directly from bottle or using a soft wipe saturated with solution. Clean and dry with soft-lint-free wipe or cloth.

- **Microscopes**
- **Optics**
- **Glass**
- **lenses**
- **Adsorption Cells**

Decon Catalog #1100

Packaging

2 oz. Squeeze Bottle
6 Bottles/pack
8 packs/case

PRECAUTIONS: Keep away from heat, sparks, or open flame. Use adequate ventilation. Keep out of reach of children. Wearing gloves and goggles is recommended when using large quantities of this product. Do not take internally. Not for ophthalmic use.

IF SWALLOWED: (if conscious) - give 1 to 2 glasses of water and call a physician. If inhaled, remove to fresh air. **FOR SKIN CONTACT:** flush with water. **FOR EYE CONTACT:** flush with water or eyewash for 15 minutes and seek medical attention.

STORAGE: Keep contents in a clean and dry environment.

www.deconlabs.com



Rev. 12/2022