



Peroxigen ST

Sterile Hydrogen Peroxide Solution



Peroxigen ST is a sterile 6% (v/v) Hydrogen Peroxide solution made with Water-for-Injection (WFI). This ready-to-use solution can be used to clean surfaces and equipment found in cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, hydrogen peroxide solution is required.

Peroxigen ST is filtered to 0.2 microns, bottled and double-bagged within a ISO CLASS 5 cleanroom. It is gamma-irradiated to a sterility assurance level (SAL) of 10^{-6} . Each lot is tested for sterility via USP 14-day sterility test. A QR code can be found on each case of **Peroxigen ST** which is linked to a copy of the lot specific document, which details QC, irradiation, and sterility for that particular lot.

INGREDIENTS

6% (v/v) Hydrogen Peroxide

94% (v/v) Water

Precautionary Statements

Causes serious eye irritation. Wear eye protection/face protection.

If in eyes: Rinse cautiously with water for several minutes.

PACKAGING

12 x 16oz.

Trigger Spray Bottle

DECON CAT NO. 9716

SHELF LIFE
2 YEARS

Rev. 12/2022