



# **Technical Data Sheet**

# -BRIEF DESCRIPTION -

Peroxigen is a ready-to-use 6% Hydrogen Peroxide solution in a convenient trigger spray bottle. Used to clean hard, non-porous surfaces in labs and production areas.

Safe to use, does not leave residues and breaks down into water and oxygen. Low-toxicity, non-flammable, and zero VOC profile helps to assist in VOC compliance.

# - INGREDIENTS -

6% (v/v) Hydrogen Peroxide 94% (v/v) Water

## - SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date is on the product.

# Plegraphrouse, (Injury Processor Law Of Action Control Control

16oz. Trigger Spray Bottle

# - PRECAUTIONARY STATEMENT -

Causes serious eye irritation.
Wear eye protection/face protection

If in eyes: Rinse cautiously with water for several minutes.

## - PACKAGING -

16 oz. (473mL)Trigger Spray / 12 per case

# - CATALOG NUMBERS -

ManufacturerItemCat. No.Decon Labs, Inc. 16 oz. Bottle9816

For Additional Technical Information call Decon Customer Service (800) 332-6647, or See our website at www.deconlabs.com

**®**Decon Labs, Inc. 06/24/2016